

# GeneProof®



## Immunocompromised / Transplanted Panel

**The widest portfolio of PCR Kits for precise, on-time and complex diagnostics.**



- + GeneProof **Cytomegalovirus (CMV) PCR Kit**
- + GeneProof **Epstein-Barr Virus (EBV) PCR Kit**
- + GeneProof **BK/JC Virus (BK/JC) PCR Kit**
- + GeneProof **BK Virus (BKV) PCR Kit**
- + GeneProof **JC Virus (JCV) PCR Kit**
- + GeneProof **Adenovirus PCR Kit**
- + GeneProof **Aspergillus PCR Kit**
- + GeneProof **Herpes Simplex Virus (HSV-1/2) PCR Kit**
- + GeneProof **Human Herpesvirus 6/7 (HHV-6/7) PCR Kit**
- + GeneProof **Human Herpesvirus 8 (HHV-8) PCR Kit**
- + GeneProof **Parvovirus B19 PCR Kit**
- + GeneProof **Varicella-Zoster Virus (VZV) PCR Kit**
- + GeneProof **Universal Internal Control**

**Combine examinations according to your laboratory needs!**

**one sample  
one extraction  
one run**



**many results**

### TECHNOLOGICALLY ADVANCED PRODUCTS

#### Easy-to-Use concept

- Universal PCR profile and Universal Internal Control simplifies laboratory workflow
- No additional pipetting of Ready-to-Use Master Mix necessary
- Many results in a single run

#### Contamination prevention

- Master Mix contains Uracil-DNA glycosylase (UNG) and dUTPs eliminating carryover contamination

#### False negative results control

- Unique solution of competitive internal control

### EXCELLENT QUALITY

- Quality proven by EQA panel results

### FULL TRACEABILITY TO W.H.O. STANDARDS



- Exact and repeatable quantification of viral load based on WHO Standards

### HIGH LEVEL OF CUSTOMER SUPPORT

- Service guaranteed within 24 hours
- Expert on the line
- Result consultation
- Laboratory workflow analysis
- Training program



**COMPATIBLE WITH A WIDE RANGE OF REAL-TIME PCR DEVICES**



**CERTIFIED DIAGNOSTIC TESTS**

[www.geneproof.com](http://www.geneproof.com)

GUARANTEED CONTROL OVER THE COMPLETE PROCESS OF DEVELOPMENT, MANUFACTURING AND DISTRIBUTION OF ALL THE OFFERED GENEPROOF PRODUCTS

# GeneProof®




## Immunocompromised / Transplanted Panel

### Complete Molecular Diagnostics Solution for Immunocompromised and Transplanted Patients

Combine GeneProof PCR Kits and make examinations according  
to your laboratory needs.

**Just multiplex!**



**Because only the complete  
solution makes sense.**

[www.geneproof.com](http://www.geneproof.com)

GeneProof a.s.

Vídeňská 101/119 / Dolní Heršpice / 619 00 Brno / Czech Republic  
+420 543 211 679 PHONE / [info@geneproof.com](mailto:info@geneproof.com) E-MAIL