

**MEM Eagle (2x) with
Earle's Salt, w/o Phenol
Red, w Sodium
Bicarbonate**

Code: *ECB2069*

Inorganic Salts	Mg/L
Calcium chloride x 2H ₂ O	530.00
Magnesium sulfate anhydrous	195.44
Potassium chloride	800.00
Sodium chloride	13,600.00
Sodium dihydrogen phosphate anhyd.	244.00

Other Components	Mg/L
D(+)-Glucose	2,000.00
Phenol red	0.00
NaHCO ₃	4,400.00

Amino Acids	Mg/L
L-Arginine x HCl	252.00
L-Cystine x 2HCl	62.60
L-Glutamine	0.00
L-Histidine x HCl x H ₂ O	84.00
L-Isoleucine	104.00
L-Leucine	104.00
L-Lysine x HCl	145.00
L-Methionine	30.00
L-Phenylalanine	64.00
L-Threonine	96.00
L-Tryptophan	20.00
L-Tyrosine disodium salt x 2H ₂ O	103.80
L-Valine	92.00

Vitamins	Mg/L
D-Calcium pantothenate	2.00
Choline chloride	2.00
Folic acid	2.00
myo-Inositol	4.00
Nicotinamide	2.00
Pyridoxal x HCl	2.00
Riboflavin	0.20
Thiamine x HCl	2.00