



EuroClone®

serving science through innovation

Dulbecco's Modified Eagle's Medium w/o Glucose, w/o L-Glutamine, w/o Sodium pyruvate, w/o Phenol Red

**Code:
ECM0062L**

Amino acids :	
L-Arginine x HCL	84,00 mg/L
L-Cystine x2HCL	62,58 mg/L
Glycine	30,00 mg/L
Histidine xHCL xH ₂ O	42,00 mg/L
L-Isoleucine	104,80 mg/L
L-Leucine	104,80 mg/L
L-Lysine x HCL	146,20 mg/L
L-Methionine	30,00 mg/L
L-Phenylalanine	66,00 mg/L
L-Serine	42,00 mg/L
L-Threonine	95,20 mg/L
L-Tryptophan	16,00 mg/L
L-Tyrosine x Na	103,79 mg/L
L-Valine	93,60 mg/L
Inorganic salt :	
Calcium chloride anhydrous	200,00 mg/L
Iron (III) nitrate x 9H ₂ O	0,10 mg/L
Magnesium sulfate anhydrous	97,66 mg/L
Potassium chloride	400,00 mg/L
Sodium chloride	6400,00 mg/L
Sodium dihydrogen phosphate anhydrous	108,69 mg/L
NaHCO ₃	3700,00 mg/L

Vitamins :	
Choline chloride	4,00 mg/L
D-Ca pantothenate	4,00mg/L
Folic acid	4,00 mg/L
Myo-Inositol	7,00 mg/L
Nicotinamide	4,00 mg/L
Pyridoxal x HCL	4,00 mg/L
Riboflavin	0,40 mg/L
Thiamine x HCL	4,00 mg/L